

# RNF156 antibody

Catalog Number: 114744

## Product name

RNF156 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RNF156 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HEK-293 cells, HeLa cells, HepG2 cells. Positive IHC detected in human lung cancer tissue. Observed molecular weight by Western-blot: 66-72 kDa, 55 kDa

## Preparation

This antibody was obtained by immunization of RNF156 recombinant protein (Accession Number: BC014199). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114744(MGRN1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:114744(MGRN1 antibody) at dilution of 1:100 (under 10x lens)