

VGLL1 antibody



Catalog Number: 116744

Product name

VGLL1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

VGLL1 Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue. Positive WB detected in Jurkat cells, human brain tissue, mouse brain tissue, mouse spleen tissue, PC-3 cells, rat brain tissue. Observed molecular weight by Western-blot: 26-29 kDa

Preparation

This antibody was obtained by immunization of VGLL1 recombinant protein (Accession Number: NM_016267). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116744(VGLL1 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:116744(VGLL1 antibody) at dilution of 1:50 (under 10x lens)