

# KIFC3 antibody

Catalog Number: 112054

## Product name

KIFC3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

KIFC3 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, mouse kidney tissue. Positive IP detected in mouse kidney tissue. Observed molecular weight by Western blot: 80-85 kDa, 95-100 kDa

## Preparation

This antibody was obtained by immunization of KIFC3 recombinant protein (Accession Number: XM\_047434084). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:112054(KIFC3 antibody) at dilution of 1:500



IP Result of anti-KIFC3 (IP:Catalog No:112054, 4ug; Detection:Catalog No:112054 1:500) with mouse kidney tissue lysate 7200ug.