

Rhesus ENPP3 / CD203c (His Tag) recombinant protein



Catalog Number: 503054

General Information

Protein Construction

A DNA sequence encoding the rhesus ENPP3 (XP_001103528.1) (Arg46-Ile874) was expressed with a polyhistidine tag at the N-terminus.

Organism

Rhesus

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant rhesus ENPP3 consists of 845 amino acids and has a calculated molecular mass of 97.2 kDa.

Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10% glycerol, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Cynomolgus ENPP3 / CD203c Protein (His Tag)
SDS-PAGE