

NDEL1 antibody

Catalog Number: 113054

Product name

NDEL1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NDEL1 Rabbit Polyclonal antibody. Positive IP detected in HepG2 cells. Positive WB detected in human brain tissue, HepG2 cells. Observed molecular weight by Western-blot: 38-42kd

Preparation

This antibody was obtained by immunization of NDEL1 recombinant protein (Accession Number: XM_017025185). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113054(NDEL1 antibody) at dilution of 1:500



IP Result of anti-NDEL1 (IP:Catalog No:113054, 4ug; Detection:Catalog No:113054 1:300) with HepG2 cells lysate 2400ug.