

# Anti-IL-1R1/CD121a antibody



Catalog Number: 104054

## Product name

Anti-IL-1R1/CD121a antibody

## Immunogen

[Rat IL-1R1/CD121a \(His Tag\) recombinant protein](#)

## Specificity

Rat IL-1R1 / CD121a

## Antibody description

Rabbit monoclonal to IL-1R1/CD121a

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Rat IL-1R1 / CD121a (rR IL-1R1 / CD121a; NP\_037255.3; Met1-Lys352) and conjugated with FITC under optimum conditions, the unreacted FITC was removed.

## Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

## Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

FCM

## Dilutions

FCM

## Validations



IL-1R1 / CD121a Antibody (FITC), Rabbit MAb,  
Flow cytometric analysis

Flow cytometric analysis of Rat IL1R1(CD121a) expression on rat splenocytes. SD splenocytes were stained with FITC-conjugated anti-Rat IL1R1(CD121a). The histogram were derived from gated events with the forward and side light-scatter characteristics of intact cells.