

Anti-PAK3 antibody



Catalog Number: 175054

Product name

Anti-PAK3 antibody

Specificity

Human, Monkey

Antibody description

Mouse monoclonal antibody to PAK3

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*TBS (pH7.4), 1%BSA, 40%Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB

Dilutions

WB: 1:500

Validations






Fig1: Western blot analysis of PAK3 on human PAK3 recombinant protein using anti-PAK3 antibody at 1/1,000 dilution.



antibody at 1/1,000 dilution.

Fig2: Western blot analysis of PAK3 on HEK293 (1) and PAK3-hIgGFc transfected HEK293 (2) cell lysate using anti-PAK3 antibody at 1/1,000 dilution.



dilution.

Fig3: Western blot analysis of PAK3 on different cell lysate using anti-PAK3 antibody at 1/1,000 dilution.; Positive control: Line1: HeLa Line2: SK-N-SH Line3: T47D Line4: COS7 Line5: HepG2