

Anti-HMGB2 Rabbit antibody



Catalog Number: 168054

Product name

Anti-HMGB2 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

HMGB2 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human HMGB2

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP

Dilutions


WB: 1/2000-1/10000

IHC: 1/50

ICC/IF: 1/100

IP: 1/20

Validations



Western blot detection of HMGB2 in CHO-K1, Ramos cell lysates using HMGB2 Rabbit mAb(1:1000 diluted). Predicted band size:24kDa. Observed band size:24kDa.