

RAB5A-Specific antibody

Catalog Number: 114454

Product name

RAB5A-Specific antibody

WB: 1:500-1:5000

IP: 1:200-1:2000

Specificity

Human, Mouse; other species not tested.


IHC: 1:20-1:200

IF: 1:20-1:200

Antibody description

RAB5A-Specific Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue, human placenta tissue. Positive IF detected in NIH/3T3 cells. Positive WB detected in human brain tissue, mouse colon tissue, SH-SY5Y cells. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 24 kDa

Validations



Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004162). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114454(RAB5A-Specific antibody) at dilution of

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

ELISA, WB, IHC, IP, IF

Dilutions


Recommended Dilution:




Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:114454(RAB5A-Specific Antibody) at dilution of 1:200 (under 10x lens).

RAB5A-Specific antibody


Catalog Number: 114454



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:114454(RAB5A-Specific Antibody) at dilution of 1:200 (under 40x lens).



Immunofluorescent analysis of (-20oc Ethanol) fixed NIH/3T3 cells using Catalog No:114454(RAB5A-Specific Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



IP Result of anti-RAB5A-Specific (IP:Catalog No:114454, 4ug; Detection:Catalog No:114454 1:600) with mouse brain tissue lysate 4000ug.