

SNRNP200 antibody



Catalog Number: 115454

Product name

SNRNP200 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

SNRNP200 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, Jurkat cells, NIH/3T3 cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 200-245 kDa

Preparation

This antibody was obtained by immunization of SNRNP200 recombinant protein (Accession Number: BC001417). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115454(SNRNP200 antibody) at dilution of



IP Result of anti-SNRNP200 (IP:Catalog No:115454, 4ug; Detection:Catalog No:115454 1:1000) with HeLa cells lysate 4000ug.