

Phospho-TP73-Y99 pAb



Catalog Number: 166454

Product name

Rabbit IgG

Phospho-TP73-Y99 pAb

Applications

Specificity

WB

Human

Dilutions

Antibody description


WB 1:500 - 1:2000

Polyclonal antibody to Phospho-TP73-Y99

Preparation

Validations

Antigen: A phospho specific peptide corresponding to residues surrounding Y99 of human TP73



Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Western blot - Phospho-TP73-Y99 pAb

Western blot analysis of extracts from 293 cells, using Phospho-TP73-Y99 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% BSA.