

AIG1 antibody

Catalog Number: 107854

Product name

AIG1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AIG1 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Positive WB detected in mouse testis tissue, mouse heart tissue, mouse ovary tissue, PC-3 cells. Observed molecular weight by Western-blot: 30-35 kDa

Preparation

This antibody was obtained by immunization of AIG1 recombinant protein (Accession Number: NM_016108). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107854(AIG1 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:107854(AIG1 antibody) at dilution of 1:100 (under 40x lens)