

# DACH2 antibody



Catalog Number: 109854

## Product name

DACH2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

DACH2 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, HeLa cells, mouse eye tissue, mouse kidney tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 65 kDa

## Preparation

This antibody was obtained by immunization of DACH2 recombinant protein (Accession Number: BC114950). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109854(DACH2 Antibody) at dilution of 1:600