

E2F3 antibody

Catalog Number: 110154

Product name

E2F3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

E2F3 Rabbit Polyclonal antibody. Positive WB detected in Y79 cells, HEK-293 cells, HepG2 cells, mouse liver tissue. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 52 kDa

Preparation

This antibody was obtained by immunization of E2F3 recombinant protein (Accession Number: BC016847). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:110154(E2F3 Antibody) at dilution of 1:600



IP Result of anti-E2F3 (IP:Catalog No:110154, 4ug; Detection:Catalog No:110154 1:500) with mouse liver tissue lysate 4000ug.