

PPP5C antibody

Catalog Number: 114154

Product name

PPP5C antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PPP5C Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive WB detected in mouse testis tissue, human liver tissue. Observed molecular weight by Western-blot: 57 kDa

Preparation

This antibody was obtained by immunization of PPP5C recombinant protein (Accession Number: NM_006247). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114154(PPP5C antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:114154(PPP5C antibody) at dilution of 1:100 (under 10x lens)