

# Anti-dars2 antibody



Catalog Number: 175154

---

## Product name

Anti-dars2 antibody

## Specificity

Zebrafish

## Antibody description

Rabbit polyclonal antibody to dars2

## Preparation

This antigen of this antibody was recombinant protein within zebrafish dars2 aa 43-321.

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 50% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

IF

## Dilutions

IF:1:50-1:200

## Validations




Fig1:; IF staining of DARS2 in whole Zebrafish tissue (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.