

# TNNT1 antibody

Catalog Number: 116154

---

**Product name**

Polyclonal

TNNT1 antibody

**Ig Type****Specificity**

Rabbit IgG

Human,Mouse,Rat; other species not tested.

**Applications****Antibody description**

ELISA, WB

TNNT1 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, human skeletal muscle tissue. Observed molecular weight by Western-blot: 30-33kd


**Dilutions**

Recommended Dilution:

WB: 1:1000-1:10000

**Preparation**

This antibody was obtained by immunization of TNNT1 recombinant protein (Accession Number: BC022086). Purification method: Antigen affinity purified.

**Validations****Formulation**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:116154(TNNT1 antibody) at dilution of 1:1500

**Clonality**