

HSDL1 antibody



Catalog Number: 111554

Product name

HSDL1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

HSDL1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Positive IHC detected in human liver. Observed molecular weight by Western-blot: 37-40kd

Preparation

This antibody was obtained by immunization of HSDL1 recombinant protein (Accession Number: NM_031463). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111554(HSDL1 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human liver using Catalog No:111554(HSDL1 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:111554(HSDL1 antibody) at dilution of 1:100 (under 40x lens)