

# MT-ATP8 Polyclonal Antibody

Catalog Number: 163554

## Product name

MT-ATP8 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to MT-ATP8

## Preparation

Antigen: Recombinant protein of human MT-ATP8.

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IF 1:50 - 1:100

## Validations




## Western blot - MT-ATP8 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using MT-ATP8 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 2min.



## Immunofluorescence - MT-ATP8 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using MT-ATP8 antibody at dilution of 1:100. Blue: DAPI



for nuclear staining.

# MT-ATP8 Polyclonal Antibody



Catalog Number: 163554

---

Immunofluorescence - MT-ATP8 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using MT-ATP8 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.