

Pan Phospho-Tyrosine pAb



Catalog Number: 166554

Product name

Pan Phospho-Tyrosine pAb

Specificity

ALL

Antibody description

Polyclonal antibody to Pan Phospho-Tyrosine

Preparation

Antigen: A synthetic peptide corresponding to a sequence containing phosphorylated S & T.

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - Pan Phospho-Tyrosine pAb

Western blot analysis of extracts of Jurkat cells, using Pan Phospho-Tyrosine pAb at 1:1000 dilution. Jurkat cells were treated by Pervanadate (1 mM) at 37°C for 30 minutes. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% BSA. Detection: ECL Basic Kit. Exposure time: 1s.