

# BUD13 antibody



Catalog Number: 108554

## Product name

BUD13 antibody

## Specificity

Human; other species not tested.

## Antibody description

BUD13 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in HepG2 cells. Observed molecular weight by Western-blot: 70-75 kDa

## Preparation

This antibody was obtained by immunization of BUD13 recombinant protein (Accession Number: NM\_032725). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:50-1:500

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:108554(BUD13 Antibody) at dilution of 1:600



Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog No:108554(BUD13 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)