

# Human HPRT1/HPRT (His tag) recombinant protein



Catalog Number: 518554

## General Information

### Gene Name Synonym

Hypoxanthine-Guanine Phosphoribosyltransferase, HGPRT, HGPRTase, HPRT1, HPRT

### Protein Construction

A DNA sequence encoding the Human HPRT1/HPRT Met1-Ala218 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

## Stability

Store for up to 6 months at -20°C

## Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris, 250mM NaCl, 2mM EDTA, 30% Glycerol, pH 8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

## SDS-PAGE



Human HPRT1/HPRT (His tag) recombinant protein