

Human PEX11A / Peroxin-11A (Fc Tag) recombinant protein



Catalog Number: 500954

General Information

Gene Name Synonym

cDNA FLJ10553 fis, clone NT2RP2002325, highly similar to Peroxisomal membrane protein 11A; cDNA, FLJ93834, Homo sapiens peroxisomal biogenesis factor 11A (PEX11A)

Protein Construction

A DNA sequence encoding the human PEX11A (NP_003838.1) (Arg106-Pro219) was expressed with the Fc region of mouse IgG1 at the N-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp

Molecular Mass

The recombinant human PEX11A consists 350 amino acids and predicts a molecular mass of 40 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human PEX11A / Peroxin-11A Protein (Fc Tag)