

Anti-CAMKV antibody



Catalog Number: 142954

Product name

Anti-CAMKV antibody

Immunogen

[Human CAMKV \(His & GST Tag\) recombinant protein](#)

Specificity

Human CAMKV

Antibody description

Rabbit Polyclonal to Human CAMKV

Preparation

Produced in rabbits immunized with purified, recombinant Human CAMKV. CAMKV specific IgG was purified by Human CAMKV affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P, ICC/IF, IF

Dilutions


ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Human CAMKV.


IHC-P: 0.1-2µg/mL

ICC/IF: 0.5-3µg/mL

Validations



Human CAMKV Immunofluorescence(IF) 15689



Human CAMKV Immunohistochemistry(IHC) 15688