

Anti-STK17B antibody



Catalog Number: 175954

Product name

Anti-STK17B antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to STK17B

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human stk17b. aa range:251-300

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of various cells using DRAK2 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventibiotech, MN, USA).