

ALAS2 antibody

Catalog Number: 107954

Product name

ALAS2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ALAS2 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of ALAS2 recombinant protein (Accession Number: NM_001037968). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107954(ALAS2 Antibody) at dilution of 1:700