

Harmonin antibody



Catalog Number: 111254

Product name

Harmonin antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Harmonin Rabbit Polyclonal antibody. Positive WB detected in human small intestine tissue. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 73kd

Preparation

This antibody was obtained by immunization of Harmonin recombinant protein (Accession Number: NM_001297764). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human small intestine tissue were subjected to SDS PAGE followed by western blot with Catalog No:111254(USH1C antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:111254(USH1C antibody) at dilution of 1:50 (under 10x lens)