

SNRPD1 Polyclonal Antibody



Catalog Number: 165254

Product name

SNRPD1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to SNRPD1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-119 of human SNRPD1 (NP_008869.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:100

Validations



Immunofluorescence - SNRPD1 Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using SNRPD1 antibody .