

Anti-CD95/APO-1/TNFRSF6 antibody



Catalog Number: 106254

Product name

Anti-CD95/APO-1/TNFRSF6 antibody

Immunogen

[Rat CD95/APO-1/TNFRSF6 \(His Tag\) recombinant protein](#)

Specificity

Rat CD95 / APO-1 / TNFRSF6

Antibody description

Rabbit monoclonal to CD95/APO-1/TNFRSF6

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Rat CD95 / APO-1 / TNFRSF6 (rR CD95 / APO-1 / TNFRSF6; NP_631933.2; Met1-Lys170) and conjugated with PE under optimum conditions, the unreacted PE was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



CD95 / APO-1 / TNFRSF6 Antibody (PE), Rabbit MAb, Flow Cytometry

Flow cytometric analysis of Rat FAS(CD95) expression on SD rat splenocytes. Cells were stained with PE-conjugated anti-Rat FAS(CD95). The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.