

TOR3A antibody



Catalog Number: 116254

Product name

TOR3A antibody

Specificity

human, mouse; other species not tested.

Antibody description

TOR3A Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 46kd

Preparation

This antibody was obtained by immunization of TOR3A recombinant protein (Accession Number: NM_022371). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116254(TOR3A antibody) at dilution of 1:300