

Anti-Nuf2 Rabbit antibody



Catalog Number: 167254

Product name

Anti-Nuf2 Rabbit antibody

Specificity

Human

Antibody description

Nuf2 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Nuf2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/1000-1/5000

ICC/IF: 1/100

FC: 1/20

Validations



Western blot detection of Nuf2 in K562 cell lysates using Nuf2 Rabbit pAb(1:1000 diluted). Predicted band size:54kDa. Observed band size:54kDa.