

# Anti-SAFB1 (2E8) Mouse antibody



Catalog Number: 168254

## Product name

Anti-SAFB1 (2E8) Mouse antibody

## Specificity

Human

## Antibody description

SAFB1 (2E8) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human SAFB1 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


WB, ICC/IF, ChIP

## Dilutions


WB: 1/1000

;ICC: 1/200

## Validations



Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using SAFB1 mouse mAb (dilution 1:200).



Western blot detection of SAFB1 in HeLaNE, A431, MDA-MB-468, Breast cancer, ZR75-1 and MCF7 cell lysates using SAFB1 mouse mAb (1:4000 diluted). Predicted band size: 130kDa. Observed band size: 130kDa.