

GSS Polyclonal Antibody



Catalog Number: 161654

Product name

GSS Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to GSS

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-474 of human GSS (NP_000169.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Western blot - GSS Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using GSS antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.



Immunohistochemistry - GSS Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using GSS antibody at dilution of



1:200 (40x lens).

GSS Polyclonal Antibody



Catalog Number: 161654

Immunofluorescence - GSS Polyclonal Antibody

Immunofluorescence analysis of A549 cells using GSS antibody . Blue: DAPI for nuclear staining.