

MEOX2 antibody

Catalog Number: 112654

Product name

MEOX2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MEOX2 Rabbit Polyclonal antibody. Positive IP detected in human placenta tissue. Positive WB detected in human placenta tissue. Observed molecular weight by Western-blot: 33 kDa

Preparation

This antibody was obtained by immunization of MEOX2 recombinant protein (Accession Number: NM_005924). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:112654(MEOX2 Antibody) at dilution of 1:600



IP Result of anti-MEOX2 (IP:Catalog No:112654, 4ug; Detection:Catalog No:112654 1:300) with human placenta tissue lysate 1520ug.