

Anti-KIAA0100 antibody



Catalog Number: 176654

Product name

Anti-KIAA0100 antibody

Specificity

Human

Antibody description

Mouse monoclonal antibody to KIAA0100

Preparation

This antigen of this antibody was recombinant protein within human kiaa0100 1900-2250 aa.

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, ICC, IHC-P

Dilutions

WB: 1:1,000-1:10,000

ICC: 1:100-1:500

IHC-P: 1:200

Validations



Fig1: Western blot analysis on recombinant protein using anti-KIAA0100 Mouse mAb (Cat. #




Fig2: Immunofluorescent staining of HeLa cells using anti-KIAA0100 Mouse mAb (Cat. #




Fig3: Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using anti-KIAA0100 Mouse mAb (Cat. # 176654#).