

# Anti-KIAA0100 antibody



Catalog Number: 176654

## Product name

Anti-KIAA0100 antibody

## Specificity

Human

## Antibody description

Mouse monoclonal antibody to KIAA0100

## Preparation

This antigen of this antibody was recombinant protein within human kiaa0100 1900-2250 aa.

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

WB, ICC, IHC-P

## Dilutions

WB: 1:1,000-1:10,000

ICC: 1:100-1:500

IHC-P: 1:200

## Validations



Fig1: Western blot analysis on recombinant protein using anti-KIAA0100 Mouse mAb (Cat. #

176654#).




Fig2: Immunofluorescent staining of HeLa cells using anti-KIAA0100 Mouse mAb (Cat. #

176654#).




Fig3: Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using anti-KIAA0100 Mouse mAb (Cat. # 176654#).