

Anti-TCP1 alpha/CCTA Rabbit antibody



Catalog Number: 167654

Product name

Rabbit IgG

Anti-TCP1 alpha/CCTA Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

TCP1 alpha/CCTA Rabbit monoclonal antibody

IHC: 1/50

Preparation

ICC/IF: 1/20


Antigen: A synthetic peptide of human TCP1 alpha/CCTA

IP: 1/20

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Validations



Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of TCP1 alpha/CCTA in HeLa,CHO-K1,Ramos using TCP1 alpha/CCTA Rabbit mAb(1:1000 diluted)

Clonality

Monoclonal

Ig Type