

# CXXC5 antibody



Catalog Number: 109654

## Product name

CXXC5 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CXXC5 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, SH-SY5Y cells. Observed molecular weight by Western-blot: 35-40 kDa

## Preparation

This antibody was obtained by immunization of CXXC5 recombinant protein (Accession Number: BC017439). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109654(CXXC5 Antibody) at dilution of 1:600