

NUDT10 antibody



Catalog Number: 113354

Product name

NUDT10 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

NUDT10 Rabbit Polyclonal antibody. Positive IF detected in Hela. Positive WB detected in human brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 23kd

Preparation

This antibody was obtained by immunization of NUDT10 recombinant protein (Accession Number: NM_001304963). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113354(NUDT10 antibody) at dilution of



Immunofluorescent analysis of Hela cells, using NUDT10 antibody Catalog No:113354 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).