

PYCR1 antibody

Catalog Number: 114354

Product name

PYCR1 antibody

IP: 1:500-1:5000

IF: 1:50-1:500

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PYCR1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, COLO 320 cells. Observed molecular weight by Western-blot: 33 kDa

Preparation

This antibody was obtained by immunization of PYCR1 recombinant protein (Accession Number: BC022244). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IP, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114354(PYCR1 antibody) at dilution of 1:1000



IP Result of anti-PYCR1 (IP:Catalog No:114354, 4ug; Detection:Catalog No:114354 1:1000) with HeLa cells lysate 3200ug.



Immunofluorescent analysis of (10% Formaldehyde) fixed HeLa cells using Catalog

PYCR1 antibody



Catalog Number: 114354

No:114354(PYCR1 Antibody) at dilution of 1:50

and Alexa Fluor 488-conjugated AffiniPure Goat
Anti-Rabbit IgG(H+L)