

Anti-GAPDH antibody



Catalog Number: 105354

Product name

Anti-GAPDH antibody

Specificity

Human GAPDH

Antibody description

Rabbit monoclonal to GAPDH

Preparation

This antibody was obtained from a rabbit immunized with A synthetic peptide corresponding to the center region of the human GAPDH.

Formulation

0.2 μ m filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions

IHC-P: 1-10 μ g/mL

Validations



GAPDH Loading Control Antibody, Rabbit Mab, Immunohistochemistry

Immunochemical staining of human GAPDH in human appendicitis with rabbit monoclonal antibody (5 μ g/mL, formalin-fixed paraffin embedded sections).