

ATP5G2-Specific antibody

Catalog Number: 108354

Product name

ATP5G2-Specific antibody

Specificity

Human; other species not tested.

Antibody description

ATP5G2-Specific Rabbit Polyclonal antibody.
Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005176).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:1000

Validations



IP Result of anti-ATP5G2-Specific (IP:Catalog No:108354, 3ug; Detection:Catalog No:108354 1:300) with HEK-293 cells lysate 3200ug.