

# Mouse Lgmn/Prsc1 (His tag) recombinant protein



Catalog Number: 518354

## General Information

### Gene Name Synonym

Legumain;Lgmn;Asparaginyl endopeptidase;Protease cysteine 1;Prsc1

### Protein Construction

A DNA sequence encoding the Mouse Lgmn/Prsc1 Val18-Tyr435 is expressed with a 6His tag at the C-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2  $\mu$ m filtered solution of 20mM Tris,150mM NaCl,20% Glycerol,pH8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Mouse Lgmn/Prsc1 (His tag) recombinant protein