

# NCF1 Polyclonal Antibody



Catalog Number: 160754

## Product name

NCF1 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to NCF1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-390 of human NCF1 (NP\_000256.4).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IF, IP


## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100


IP 1:50 - 1:100

## Validations



Western blot - NCF1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using NCF1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunoprecipitation - NCF1 Polyclonal Antibody

Immunoprecipitation analysis of 200ug extracts of Raji cells, using 3 ug NCF1 antibody. Western blot was performed from the immunoprecipitate using NCF1 antibody at a dilution of 1:1000.