

Vimentin antibody



Catalog Number: 116754

Product name

Polyclonal

Vimentin antibody

Ig Type

Specificity

Rabbit IgG

Human,Mouse,Rat; other species not tested.

Applications

Antibody description

ELISA, WB

Vimentin Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, mouse heart tissue. Observed molecular weight by Western-blot: 57kd

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of Vimentin recombinant protein (Accession Number: XM_006717500). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116754(VIM antibody) at dilution of 1:1000

Clonality