

Anti-phospho-Cdk2 (Thr14) Rabbit antibody



Catalog Number: 167754

Product name

Anti-phospho-Cdk2 (Thr14) Rabbit antibody

Specificity

Human, Rat

Antibody description

phospho-Cdk2 (Thr14) Rabbit monoclonal antibody

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Thr14 of human Cdk2

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP

Dilutions


WB: 1/1000

IHC: 1/50

ICC/IF: 1/50

IP: 1/20

Validations



Western blot detection of Cdk2 (phospho-Thr14) in HeLa,A549,HL-60,U2OS using Cdk2 (phospho-Thr14) Rabbit mAb(1:1000 diluted)