

# NDUFA8 antibody



Catalog Number: 113064

## Product name

NDUFA8 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

NDUFA8 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, mouse heart tissue, mouse liver tissue. Observed molecular weight by Western-blot: 19kd

## Preparation

This antibody was obtained by immunization of NDUFA8 recombinant protein (Accession Number: NM\_014222). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:113064(NDUFA8 antibody) at dilution of 1:500