

Anti-MAGEA10-MAGEA5 antibody



Catalog Number: 176064

Product name

Anti-MAGEA10-MAGEA5 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to MAGEA10-MAGEA5

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human magea5. aa range:68-117

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of 3T3 cells using MAGE-A5 Polyclonal Antibody