

APOD antibody



Catalog Number: 107064

Product name

APOD antibody

Specificity

Human; other species not tested.

Antibody description

APOD Mouse Monoclonal antibody. Positive WB detected in human plasma tissue. Observed molecular weight by Western-blot: 30-33 kDa

Preparation

This antibody was obtained by immunization of APOD recombinant protein (Accession Number: NM_001647). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:5000-1:50000

Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:107064(APOD Antibody) at dilution of 1:20000