

# Anti-Ferritin Heavy Chain Rabbit antibody



Catalog Number: 169064

## Product name

Anti-Ferritin Heavy Chain Rabbit antibody

## Specificity

Human

## Antibody description

Ferritin Heavy Chain Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human Ferritin.

## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal


## Ig Type

Rabbit IgG

## Applications

ELISA

## Validations



Standard Curve for Ferritin: Capture Antibody Mouse mAb (201067, Ferritin Light Chain Mouse mAb) to Ferritin at 2u03bcg/ml and Detector Antibody Rabbit pAb(3010667, Ferritin rabbit pAb) to Ferritin at 0.5u03bcg/ml.